


**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*“Bringing plant tissue culture to the classroom and home since 1998”*


**NOTE: All protocols are intended to be used following the methods described in the KCK Manual: “Plant Tissue Culture for the Classroom and Home” or the online workshop handout.**

## Plant Hormone/Growth Regulator Summary


HORMONE	MW	Milligrams in 1 $\mu\text{M}$	Milligrams in 2 $\mu\text{M}$	Milligrams in 4 $\mu\text{M}$	1 mg in $\mu\text{M}$
BAP	225	.225	.45	.90	4.44 $\mu\text{M}$
KINETIN	215	.215	.43	.86	4.65 $\mu\text{M}$
2iP	203	.203	.41	.82	4.82 $\mu\text{M}$
TDZ	220	.220	.44	.88	4.55 $\mu\text{M}$
Zeatin	219	.219	.44	.88	4.55 $\mu\text{M}$
IAA	175	.175	.35	.70	5.71 $\mu\text{M}$
IBA	203	.203	.40	.80	4.92 $\mu\text{M}$
NAA	186	.186	.37	.74	5.37 $\mu\text{M}$
24D	221	.221	.44	.88	4.55 $\mu\text{M}$
GA3	346	.366	.73	1.46	2.89 $\mu\text{M}$

### Plant hormones/growth regulators vary in their stability.

- ❖ Powdered forms of BAP and NAA are stored at room temperature.
- ❖ All other powders should be stored frozen.
- ❖ It is recommended that all powdered forms be shipped via priority mail.
- ❖ Solutions vary in their stability.
  - BAP, kinetin and NAA solutions appear to be stable at room temperature for less than 30 days but should be refrigerated as soon as possible.
  - IAA is not stable in light and will degrade within 3 days. Make fresh and keep refrigerated.
  - Other solutions need to be refrigerated or frozen. Shipping of these solutions should be avoided unless shipped within 3 days. These are best purchased in powder form and made in solution on site followed by refrigeration.



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## **BAP Dilutions**

The molecular weight of BAP, benzylaminopurine, is 225


### **Molarity per liter**

1 M	=	225 grams
0.1 M	=	22.5 grams
0.01 M	=	2.25 grams
0.001 M (1 mM)	=	0.225 grams (225 mg)
0.0001 M (0.1 mM)	=	0.0225 grams (22.5 mg)
0.00001 M (0.01 mM)	=	0.00225 grams (2.25 mg)
0.000001 M (1 $\mu$ M)	=	0.000225 grams (0.225 mg)


***NOTE: Therefore, a solution that is 2.22  $\mu$ M BAP is 2.22 \* 0.225 mg BAP = 0.4495 mg BAP, and this is rounded to 0.5 mg BAP.***

If the BAP stock solution = 1 mg/ml, then using 0.5 ml = 0.5 mg

<b>Desired concentration of BAP</b>	<b>Volume of BAP stock solution (1 mg/ml) to add for 1 liter medium</b>
1 $\mu$ M = 0.225 mg	1 ml contains 1 mg BAP so 0.225 ml contains 0.225 mg BAP
2 $\mu$ M = 0.45 mg	= 0.450 ml
3 $\mu$ M = 0.675 mg	= 0.675 ml
4 $\mu$ M = 0.90 mg	= 0.900 ml
4.44 $\mu$ M = 0.99 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## KINETIN Dilutions

The molecular weight of kinetin, 6-Furfurylaminopurine, is 215.21


### Molarity per liter

1 M	=	215.21 grams
0.1 M	=	21.52 grams
0.01 M	=	2.152 grams
0.001 M (1 mM)	=	0.2152 grams (215 mg)
0.0001 M (0.1 mM)	=	0.0215 grams (21.5 mg)
0.00001 M (0.01 mM)	=	0.00215 grams (2.15 mg)
0.000001 M (1 $\mu$ M)	=	0.000215 grams (0.215 mg)


**NOTE:** Therefore, a solution that is 2.22  $\mu$ M kinetin is calculated:  
 $2.22 * 0.215 \text{ mg kinetin} = 0.4773 \text{ mg kinetin}$ , and  
this is rounded to 0.5 mg kinetin.

If the kinetin stock solution = 1 mg/ml, then using 0.5 ml = 0.5 mg

Desired concentration of kinetin	Volume of kinetin stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.215 mg	1 ml contains 1 mg kinetin Therefore 0.215 ml contains 0.215 mg kinetin
2 $\mu$ M = 0.430 mg	= 0.43 ml
3 $\mu$ M = 0.645 mg	= 0.65 ml
4 $\mu$ M = 0.860 mg	= 0.86 ml
4.65 $\mu$ M = 0.9998 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## **TDZ Dilutions**

The molecular weight of TDZ, Thidiazuron [1-phenyl-3-(1, 2, 3-thiazol-5-yl) urea], is 220.25


### **Molarity per liter**

1 M	=	220.25 grams
0.1 M	=	22.03 grams
0.01 M	=	2.20 grams
0.001 M (1 mM)	=	0.220 grams (220 mg)
0.0001 M (0.1 mM)	=	0.0220 grams (22.0 mg)
0.00001 M (0.01 mM)	=	0.00220 grams (2.20 mg)
0.000001 M (1 $\mu$ M)	=	0.000220 grams (0.220 mg)


***NOTE: Therefore, a solution that is 2.22  $\mu$ M TDZ is 2.22 \* 0.220 mg TDZ = 0.4884 mg, and this is rounded to 0.5 mg TDZ.***

If the TDZ stock solution = 1 mg/ml, then using 0.5 ml = 0.5 mg

<b>Desired concentration of TDZ</b>	<b>Volume of TDZ stock solution (1 mg/ml) to add for 1 liter medium</b>
1 $\mu$ M = 0.220 mg	1 ml contains 1 mg TDZ Therefore 0.220 ml contains 0.220 mg TDZ
2 $\mu$ M = 0.44 mg	= 0.44 ml
3 $\mu$ M = 0.66 mg	= 0.66 ml
4 $\mu$ M = 0.88 mg	= 0.88 ml
4.55 $\mu$ M = 1.00 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## Zeatin Dilutions

The molecular weight of zeatin is 219.25

### Molarity per liter

1 M	=	219.25 grams
0.1 M	=	21.93 grams
0.01 M	=	2.19 grams
0.001 M (1 mM)	=	0.219 grams (219 mg)
0.0001 M (0.1 mM)	=	0.0219 grams (21.9 mg)
0.00001 M (0.01 mM)	=	0.0022 grams (2.19 mg)
0.000001 M (1 $\mu$ M)	=	0.000219 grams (0.219 mg)


**NOTE:** Therefore, a solution that is 2.22  $\mu$ M zeatin is  $2.22 * 0.219 \text{ mg} = 0.4862 \text{ mg}$ , and this is rounded to 0.5 mg zeatin.

If the zeatin stock solution = 1 mg/ml, then using 0.4 ml = 0.4 mg

Desired concentration of zeatin	Volume of zeatin stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.219 mg	1 ml contains 1 mg zeatin Therefore 0.219 ml contains 0.219 mg zeatin
2 $\mu$ M = 0.43 mg	= 0.43 ml
3 $\mu$ M = 0.65 mg	= 0.65 ml
4 $\mu$ M = 0.86 mg	= 0.86 ml
4.55 $\mu$ M = 0.999 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## 2iP

The molecular weight of 2iP, 6-(y,y-DIMETHYLALLYLAMINO)PURINE, is 203.3


### Molarity per liter

1 M	=	203.3 grams
0.1 M	=	20.33 grams
0.01 M	=	2.03 grams
0.001 M (1 mM)	=	0.203 grams (203 mg)
0.0001 M (0.1 mM)	=	0.0203 grams (20.3 mg)
0.00001 M (0.01 mM)	=	0.00203 grams (2.03 mg)
0.000001 M (1 $\mu$ M)	=	0.000203 grams (0.203 mg)


**NOTE:** Therefore, a solution that is 2.22  $\mu$ M 2iP is  $2.22 * 0.203 \text{ mg 2iP} = 0.4507 \text{ mg}$ , and this is rounded to 0.5 mg 2iP

If the 2iP stock solution = 1 mg/ml, then using 0.5 ml = 0.5 mg

Desired concentration of 2iP	Volume of 2iP stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.203 mg	1 ml contains 1 mg 2iP Therefore 0.203 ml contains 0.203 mg 2iP
2 $\mu$ M = 0.41 mg	= 0.41 ml
3 $\mu$ M = 0.61 mg	= 0.61 ml
4 $\mu$ M = 0.82 mg	= 0.82 ml
4.821 $\mu$ M = 0.999 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## IAA Dilutions

The molecular weight of IAA, indole-3-acetic acid, is 175.19


### Molarity per liter

1 M =	175.19 grams
0.1 M =	17.52 grams
0.01 M =	1.75 grams
.001 M (1 mM) =	0.175 grams (175 mg)
0001 M (0.1 mM) =	0.0175 grams (17.5 mg)
0.00001 M (0.01 mM) =	0.00175 grams (1.75 mg)
0.000001 M (1 $\mu$ M) =	0.000175 grams (0.175 mg)


**NOTE: Therefore, a solution that is 2.22  $\mu$ M IAA is 2.22 \* 0.175 mg IAA = 0.3885 mg, and this is rounded to 0.4 mg IAA**

If the IAA stock solution = 1 mg/ml, then using 0.4 ml = 0.4 mg

Desired concentration of IAA	Volume of IAA stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.175 mg	1 ml contains 1 mg IAA Therefore 0.175 ml contains 0.175 mg IAA
2 $\mu$ M = 0.35 mg	= 0.35 ml
3 $\mu$ M = 0.525 mg	= 0.53 ml
4 $\mu$ M = 0.70 mg	= 0.70 ml
5.71 $\mu$ M = 0.999 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## IBA Dilutions

The molecular weight of IBA, indole-3-butyric acid is 203.24


### Molarity per liter

1 M	=	203.24 grams
0.1 M	=	20.34 grams
0.01 M	=	2.03 grams
0.001 M (1 mM)	=	0.203 grams (203 mg)
0.0001 M (0.1 mM)	=	0.0203 grams (20.3 mg)
0.00001 M (0.01 mM)	=	0.00203 grams (2.03 mg)
0.000001 M (1 $\mu$ M)	=	0.000203 grams (0.203 mg)


**NOTE:** Therefore, a solution that is 2.22  $\mu$ M IBA is  $2.22 * 0.203$  mg IBA = 0.4507 mg, and this is rounded to 0.5 mg IBA.

If the IBA stock solution = 1 mg/ml, then using 0.4 ml = 0.4 mg

Desired concentration of IBA	Volume of IBA stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.203 mg	1 ml contains 1 mg IBA Therefore 0.203 ml contains 0.203 mg IBA
2 $\mu$ M = 0.41 mg	= 0.41 ml
3 $\mu$ M = 0.61 mg	= 0.61 ml
4 $\mu$ M = 0.82 mg	= 0.82 ml
4.92 $\mu$ M = 0.999 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## NAA Dilutions


The molecular weight of NAA, naphthalene acidic acid, is 186.2

Molarity per liter	
1 M	= 186.2 grams
0.1 M	= 18.62 grams
0.01 M	= 1.86 grams
0.001 M (1 mM)	= 0.186 grams (186 mg)
0.0001 M (0.1 mM)	= 0.0186 grams (18.6 mg)
0.00001 M (0.01 mM)	= 0.00186 grams (1.86 mg)
0.000001 M (1 $\mu$ M)	= 0.000186 grams (0.186 mg)


**NOTE:** Therefore, a solution that is 2.22  $\mu$ M NAA is:  $2.22 * 0.186 \text{ mg NAA} = 0.4129 \text{ mg}$ , and this is rounded to 0.4 mg NAA.

If the NAA stock solution = 1 mg/ml, then using 0.4 ml = 0.4 mg

Desired concentration of NAA	Volume of NAA stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.186 mg	1 ml contains 1 mg NAA Therefore 0.186 ml contains 0.186 mg NAA
2 $\mu$ M = 0.370 mg	= 0.37 ml
3 $\mu$ M = 0.558 mg	= 0.56 ml
4 $\mu$ M = 0.740 mg	= 0.74 ml
5.37 $\mu$ M = 0.999 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## **2,4-D Dilutions**

The molecular weight of 24D, 2,4-dichlorophenoxyacetic acid, is 221.04


### **Molarity per liter**

1 M	=	221.04 grams
0.1 M	=	22.1 grams
0.01 M	=	2.21 grams
0.001 M (1 mM)	=	0.221 grams (221 mg)
0.0001 M (0.1 mM)	=	0.0221 grams (22.1 mg)
0.00001 M (0.01 mM)	=	0.00221 grams (2.21 mg)
0.000001 M (1 $\mu$ M)	=	0.000221 grams (0.221 mg)


***NOTE: Therefore, a solution that is 2.22  $\mu$ M 24D is 2.22 \* 0.221 mg TDZ = 0.4906 mg, and this is rounded to 0.5 mg 24D.***

If the 24D stock solution = 1 mg/ml, then using 0.5 ml = 0.5 mg

<b>Desired concentration of 24D</b>	<b>Volume of TDZ stock solution (1 mg/ml) to add for 1 liter medium</b>
1 $\mu$ M = 0.221 mg	1 ml contains 1 mg 24D Therefore 0.220 ml contains 0.220 mg 24D
2 $\mu$ M = 0.44 mg	= 0.44 ml
3 $\mu$ M = 0.66 mg	= 0.66 ml
4 $\mu$ M = 0.88 mg	= 0.88 ml
4.55 $\mu$ M = 1.00 mg	= 1.000 ml



**Kitchen Culture Kits Inc.**  
Protocols for the Hobbyist  
[www.kitchenculturekit.com](http://www.kitchenculturekit.com)  
[carolstiff@kitchenculturekit.com](mailto:carolstiff@kitchenculturekit.com)



*"Bringing plant tissue culture to the classroom and home since 1998"*

## **GA3 Dilutions**

The molecular weight of GA3, gibberellic acid 3, is 346.4

1 M	=	346.4 gram
0.1 M	=	34.6 grams
0.01 M	=	3.46 grams
0.001 M (1 mM)	=	0.346 grams (346 mg)
0.0001 M (0.1 mM)	=	0.0346 grams (34.6 mg)
0.00001 M (0.01 mM)	=	0.00346 grams (3.46 mg)
0.000001 M (1 $\mu$ M)	=	0.000346 grams (0.346 mg)

**NOTE:** Therefore, a solution that is 2.22  $\mu$ M GA3 is:  $2.22 * 0.386 \text{ mg GA3} = 0.8492 \text{ mg}$ , and this is rounded to 0.9 mg GA3.

If the GA3 stock solution = 1 mg/ml, then using 0.9 ml = 0.9 mg

**\*\*NOTE: GA3 stock solution is often made at 13 mg/ml\*\***

Desired concentration of GA3	Volume of GA3 stock solution (1 mg/ml) to add for 1 liter medium
1 $\mu$ M = 0.346 mg	1 ml contains 1 mg GA3 Therefore 0.346 ml contains 0.346 mg GA3
2 $\mu$ M = 0.692 mg	= 0.69 ml
3 $\mu$ M = 1.038 mg	= 1.04 ml
4 $\mu$ M = 1.384 mg	= 1.38 ml
2.89 $\mu$ M = 0.999 mg	= 1.000 ml